

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Uniprot ID | Q9NPF7 |
| Common Name | ABBV 066;ABBV066;ABBV 066;ABBV 66;ABBV-66;ABBV66;BI-655066;BI655066;BI655066;SKYRIZ |
| Conjugate | Biotinylated |
| Synonyms | P19;SGRF;IL-23;IL-23A;IL23P19 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | IL23A |
| Description | Biotinylated Anti-IL23A(risankizumab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

