

PRODUCT INFORMATION

Uniprot ID Q13421

MORAb-009 **Common Name** Conjugate Biotinylated

MSLN; Mesothelin; MPF **Synonyms**

Applications ELISA; Flow Cyt

Recommended

ELISA 1:5000-10000; Flow Cyt 1:100 **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species Chimeric lgG1 IgG type Reactivity Human **MSLN Target**

Biotinylated Anti-Mesothelin (amatuximab biosimilar) mAb **Description**

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Email: info@dimabio.com Website: www.dimabio.com

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

