

## **PRODUCT INFORMATION**

**Uniprot ID** Q15116

**Common Name** HLX-10; HLX10; HLX 10

Conjugate Biotinylated

PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1 **Synonyms** 

**Applications** N/A

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** N/A IgG4 IgG type Reactivity Human PD-1 **Target** 

Biotinylated Anti-PD-1(serplulimab biosimilar) **Description** 

mAb

**Delivery** 6-8 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

