

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Uniprot ID | Q15116 |
| Common Name | AB 122;AB-122;AB122;GLS 010;GLS010;GLS010;GLS10;GLS-10;GSL 10;GS-0122;GS0122;GS 0122;GS-122;GS122;GS 122;WBP 3055;WBP3055;WBP-3055;YUTUO |
| Conjugate | Biotinylated |
| Synonyms | PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | PD-1 |
| Description | Biotinylated Anti-PD-1(zimberelimab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

