

PRODUCT INFORMATION

Uniprot ID Q9NZQ7

HTI-1088;HTI1088;HTI-1316;HTI1316;SHR **Common Name**

1316;SHR1316;SHR-1316

Conjugate Biotinylated

B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1 **Synonyms**

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before Formulation &

lyophilization. Please see Certificate of Analysis for Reconstitution

specific instructions of reconstitution.

Host Species N/A IgG type IgG4 Reactivity Human PD-L1 **Target**

Biotinylated Anti-PD-L1(adebrelimab biosimilar) **Description**

mAb

Delivery 6-8 weeks

Email: info@dimabio.com Website: www.dimabio.com

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid Storage & Shipping

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

