PRODUCT INFORMATION

Uniprot ID 000206

EB 05;EB05;EB-05;EB

5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB 7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101 **Common Name**

Conjugate Biotinylated

TOLL;CD284;TLR-4;ARMD10 **Synonyms**

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before Formulation & Reconstitution lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human TLR4 **Target**

Biotinylated Anti-TLR4(paridiprubart biosimilar) **Description**

mAb

6-8 weeks **Delivery**

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

