

PRODUCT INFORMATION

Uniprot ID	O00206
Common Name	EB 05;EB05;EB-05;EB 5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB 7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101
Conjugate	Biotinylated
Synonyms	TOLL;CD284;TLR-4;ARMD10
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	Human IgG1(K214R,N325S,L328F) - Kappa
Reactivity	Human
Target	TLR4
Description	Biotinylated Anti-TLR4(paridiprubart biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

