

PRODUCT INFORMATION

Target	OSMR
Synonyms	OSMRB, PLCA1, IL-31RB, OSMRbeta, IL-31R-beta
Description	Recombinant Cynomolgus OSMR protein with C-terminal 10×His tag
Delivery	In Stock
Uniprot ID	A0A2K5UFW5
Expression Host	HEK293
Tag	C-10×His tag
Molecular Characterization	OSMR(Glu28-Ser737) 10×His tag
Molecular Weight	The protein has a predicted molecular mass of 82.3 kDa after removal of the signal peptide.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a member of the type I cytokine receptor family. The encoded protein heterodimerizes with interleukin 6 signal transducer to form the type II oncostatin M receptor and with interleukin 31 receptor A to form the interleukin 31 receptor, and thus transduces oncostatin M and interleukin 31 induced signaling events. Mutations in this gene have been associated with familial primary localized cutaneous amyloidosis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]
Usage	Research use only





Figure 1. Cynomolgus OSMR Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

