

PRODUCT INFORMATION

Target	Beta-NGF
Synonyms	Beta-Nerve Growth Factor;Beta-NGF;NGF;NGFB;-NGF
Description	Recombinant Human Beta-Nerve Growth Factor is produced by our Mammalian expression system and the target gene encoding Ser122-Arg239 is expressed.
Delivery	In Stock
Uniprot ID	P01138
Expression Host	HEK293
Tag	
Molecular Characterization	Not available
Molecular Weight	13.3 KDa
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Formulation & Reconstitution	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 250mM NaCl, pH 7.0.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Human β-Nerve Growth Factor (β-NGF) was initially isolated in the mouse submandibular gland. It is composed of three non-covalently linked subunits α, β, and γ; it exhibits all the biological activities ascribed to NGF. It is structurally related to BDNF, NT-3 and NT-4 and belongs to the cysteine-knot family of growth factors that assume stable dimeric structures. B-NGF is a neurotrophic factor that signals through its receptor β-NGF, and plays a crucial role in the development and preservation of the sensory and sympathetic nervous systems. B-NGF also acts as a growth and differentiation factor for B lymphocytes and enhances B-cell survival. These results suggest that β-NGF is a pleiotropic cytokine, which in addition to its neurotropic activities may have an important role in the regulation of the immune system. Human β-NGF shares 90% sequence similarity with mouse protein and shows cross-species reactivity.
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Greater than 95% as determined by reducing SDS-PAGE.

