

PRODUCT INFORMATION

Target	EGFR
Synonyms	ERBB, ERRP, HER1, mENA, ERBB1, NNCIS, PIG61, NISBD2
Description	Recombinant human EGFR(313-445) Protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P00533
Expression Host	HEK293
Tag	C-Human Fc tag
Molecular Characterization	EGFR(Ala313-Ile445) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 41.0 kDa after removal of the signal peptide.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor, thus inducing receptor dimerization and tyrosine autophosphorylation leading to cell proliferation. Mutations in this gene are associated with lung cancer. EGFR is a component of the cytokine storm which contributes to a severe form of Coronavirus Disease 2019 (COVID-19) resulting from infection with severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2). [provided by RefSeq, Jul 2020]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human EGFR(313-445) Protein, hFc Tag on SDS-PAGE under reducing condition.

Human EGFR(313-445), hFc tagged protein ELISA

0.2 μ g of Human EGFR(313-445), hFc tagged protein per well




Figure 2. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human EGFR(313-445) Protein, hFc Tag (PME101867) can bind Anti-EGFR(depatuxizumab biosimilar) mAb (BME100252) in a linear range of 200-1000 ng/mL. In order to specifically detect BME100252, mouse anti-human Fab-specific antibody was used as detection antibody.

