

PRODUCT INFORMATION

Target	Her2
Synonyms	ERBB2;CD340;HER-2/neu;HER2;MLN19;NEU;NGL;TKR1
Description	Recombinant Human Her2 with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P04626
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	Her2(Thr23-Thr652) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 95.5 kDa after removal of the signal peptide. The apparent molecular mass of Her2-hFc is approximately 100-130 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a member of the epidermal growth factor (EGF) receptor family of receptor tyrosine kinases. This protein has no ligand binding domain of its own and therefore cannot bind growth factors. However, it does bind tightly to other ligand-bound EGF receptor family members to form a heterodimer, stabilizing ligand binding and enhancing kinase-mediated activation of downstream signalling pathways, such as those involving mitogen-activated protein kinase and phosphatidylinositol-3 kinase. Allelic variations at amino acid positions 654 and 655 of isoform a (positions 624 and 625 of isoform b) have been reported, with the most common allele, Ile654/Ile655, shown here. Amplification and/or overexpression of this gene has been reported in numerous cancers, including breast and ovarian tumors. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized.
Usage	Research use only
Conjugate	Unconjugated






Figure 1 . Human Her2 Protein, hFc Tag on SDS-PAGE under reducing condition.

**Human Her2, hFc tagged protein ELISA**  
0.1  $\mu$ g of Human Her2, hFc tagged protein per well




Figure 2. ELISA plate pre-coated by 1  $\mu$ g/mL (100  $\mu$ L/well) Human Her2 Protein, hFc Tag (PME100665) can bind Anti-Her2 (trastuzumab biosimilar) mAb (BME100048) in a linear range of 3.81-1730.10 ng/mL.



**Human Her2, hFc Tagged protein ELISA**  
0.2 µg of Human Her2, hFc tagged protein per well




Figure 3. ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human Her2 Protein, hFc Tag (PME100665) can bind Anti-HER2 antibody(DMC369), IgG1 Chimeric mAb in a linear range of 3.20-80 ng/mL.

