

**PRODUCT INFORMATION**

Target	IFN
Synonyms	Interferon Gamma;IFN-Gamma;Immune Interferon;IFNG
Description	Recombinant Human Interferon Gamma is produced by our Mammalian expression system and the target gene encoding Gln24-Gln166 is expressed.
Delivery	In Stock
Uniprot ID	P01579
Expression Host	HEK293
Molecular Characterization	Not available
Molecular Weight	16.8 KDa
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Formulation & Reconstitution	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 4% Mannitol, 2% Sucrose, 0.02% Tween80, pH7.4.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	IFN γ is the major interferon produced by mitogenically or antigenically stimulated lymphocytes. It is structurally different from type I interferon and its major activity is immunoregulation. It has been implicated in the expression of class II histocompatibility antigens in cells that do not normally produce them, leading to autoimmune disease. Interferon gamma is produced mainly by T-cells and natural killer cells activated by antigens, mitogens, or alloantigens. It is produced by lymphocytes expressing the surface antigens CD4 and CD8. IFN γ synthesis is induced by IL-2, FGF-basic, and EGF.
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Greater than 95% as determined by reducing SDS-PAGE.

DIMABIO CONFIDENTIAL

