

PRODUCT INFORMATION

Tag	C-Flag Tag
Target	KISSR
Synonyms	AXOR12; CPPB1; GPR54; HH8; HOT7T175; KISS-1R
Description	Human KISSR full length protein-synthetic nanodisc
Delivery	In Stock
Uniprot ID	Q969F8
Expression Host	HEK293
Protein Families	Druggable Genome, GPCR, Transmembrane
Protein Pathways	Neuroactive ligand-receptor interaction
Molecular Weight	The human full length KISSR protein has a MW of 42.6 kDa Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.
Formulation & Reconstitution	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Storage&Shipping	
Background	The protein is a galanin-like G protein-coupled receptor that binds metastatin, a peptide encoded by the metastasis suppressor gene KISS1. The tissue distribution of the expressed gene suggests that it is involved in the regulation of endocrine function, and this is supported by the finding that this gene appears to play a role in the onset of puberty. Mutations in this gene have been associated with hypogonadotropic hypogonadism and central precocious puberty.
Usage	Research use only
Conjugate	Unconjugated



ELISA assay to evaluate KISSR-Nanodisc 0.2 μ g Human KISSR-Nanodisc per well




Figure1. Elisa plates were pre-coated with Flag Tag KISSR-Nanodisc (0.2 μ g/per well). Serial diluted anti-Flag monoclonal antibody solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-Flag monoclonal antibody binding with KISSR-Nanodisc is 6.559ng/ml.

kDa M R

250
130
100
70
55
35
25
15
10



Figure2. WB analysis of Human KISSR-Nanodisc with anti-Flag monoclonal antibody at 1/5000 dilution, followed by Goat Anti-Rabbit IgG HRP at 1/5000 dilution

