

PRODUCT INFORMATION

Target	MXRA8
Synonyms	ASP3
Description	Recombinant Human MXRA8 Protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	Q9BRK3
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	MXRA8(Val20-Gln340) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 36.3 kDa after removal of the signal peptide. The apparent molecular mass of MXRA8-His is approximately 35-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Predicted to be involved in establishment of glial blood-brain barrier. Located in extracellular exosome. [provided by Alliance of Genome Resources, Apr 2022]
Usage	Research use only
Conjugate	Unconjugated






Figure 1.Human MXRA8 Protein, His Tag on SDS-PAGE under reducing condition.

