

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	NK2R
Synonyms	NKNAR; SKR; TAC2R
Description	Human NK2R-Strep full length protein-synthetic nanodisc
Delivery	In Stock
Uniprot ID	P21452
Expression Host	HEK293
Protein Families	Druggable Genome, GPCR, Transmembrane
Protein Pathways	Calcium signaling pathway, Neuroactive ligand-receptor interaction
Molecular Weight	The human full length NK2R-Strep protein has a MW of 44.4 kDa Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.
Formulation & Reconstitution	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Storage&Shipping	
Background	This protein belongs to a family of proteins that function as receptors for tachykinins. Receptor affinities are specified by variations in the 5'-end of the sequence. The receptors belonging to this family are characterized by interactions with G proteins and 7 hydrophobic transmembrane regions.
Usage	Research use only
Conjugate	Unconjugated



ELISA assay to evaluate NK2R-Strep-Nanodisc 0.2 μ g Human NK2R-Strep-Nanodisc per well




Figure 1. Elisa plates were pre-coated with C-Flag&Strep Tag NK2R-Strep-Nanodisc (0.2 μ g/per well). Serial diluted anti-Flag monoclonal antibody solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-Flag monoclonal antibody binding with NK2R-Strep-nanodisc is 2.571ng/ml.

kDa M R



Figure 2. Human NK2R-Strep-Nanodisc, C-Flag&Strep Tag on SDS-PAGE

