

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	NMUR1
Synonyms	FM-3, FM3, GPC-R, GPR66, NMU1R
Description	Human NMUR1-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q9HB89
Expression Host	HEK293
Protein Families	GPCR,Transmembrane,Druggable Genome,
Protein Pathways	GPCRDB Class A Rhodopsin-like,
Molecular Weight	The human full length NMUR1-Strep protein has a MW of 47.4 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Receptor for the neuromedin-U and neuromedin-S neuropeptides.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

