

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	OPSD
Synonyms	CSNBAD1, OPN2, RP4
Description	Human OPSD-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	P08100
Expression Host	HEK293
Protein Families	Transmembrane,Druggable Genome,
Protein Pathways	GPCRDB Class A Rhodopsin-like,Integrin-mediated cell adhesion KEGG,G-Protein Coupled Receptors Signaling Pathway,
Molecular Weight	The human full length OPSD-Strep protein has a MW of 38.9 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	The protein encoded by this gene is found in rod cells in the back of the eye and is essential for vision in low-light conditions. The encoded protein binds to 11-cis retinal and is activated when light hits the retinal molecule. Defects in this gene are a cause of congenital stationary night blindness. [provided by RefSeq, Aug 2017]
Usage	Research use only
Conjugate	Unconjugated

