

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	PK2L2
Synonyms	TRPP5
Description	Human PK2L2-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q9NZM6
Expression Host	HEK293
Protein Families	Ion Channels: Transient receptor potential
Protein Pathways	N/A
Molecular Weight	The human full length PK2L2-Strep protein has a MW of 73.8 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	May function as a subunit of a cation channel and play a role in fertilization.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

