

PRODUCT INFORMATION

Target	B4GALT1, B7H6, CD166, CD63, CXCR4, EPHA2, FOLR1, IGF-1R, IGFBP7, LGALS1, MSLN
Description	The embedded and sliced tumor sections (4μm section) from M-NSG mice implanted with Hela cells
Host Cells	Hela
Applications	IHC
Host mice	M-NSG
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping	Room temperature
Usage	Research use only

