

PRODUCT INFORMATION

Target	CD147
Description	The embedded and sliced tumor sections (4μm section) from M-NSG mice implanted with LnCap cells
Host Cells	LnCap
Applications	IHC
Host mice	M-NSG
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping	Room temperature
Usage	Research use only




Figure 1. IHC analysis of M-NSG LnCap DiSliceX™ SlideSet section (Cat. No. SLI100015) with Anti-CD99 antibody(DM203)Rabbit mAb (Cat. No. DME100203) at a concentration of 10 μg/ml. The inset image shows the absence of staining in secondary antibody-only control.

Address: Wuhan Institute of
Biotechnology B7, Biolake No.666
Gaoxin Road, Wuhan, Hubei, China
Telephone: +1 2409940618(USA)
/+86-18062749453(China)
/+86-400-006-0995(China)

Email: info@dimabio.com
Website: www.dimabio.com



DIMABIO CONFIDENTIAL

