

PRODUCT INFORMATION

Target	ALPP, CD63, LGALS1, YAP1
Description	The embedded and sliced tumor sections (4µm section) from M-NSG mice implanted with SiHa cells
Host Cells	SiHa
Applications	IHC
Host mice	M-NSG
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage&Shipping	Room temperature
Usage	Research use only

