

**PRODUCT INFORMATION**

<b>Target</b>	TENM4
<b>Synonyms</b>	Doc4;ELM2;l(7)-3Rn;l7Rn3;mKIAA1302;Odz4;R75022;Ten-m4
<b>Description</b>	Recombinant mouse TENM4 protein with N-terminal 6×His tag
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	Q3UHK6
<b>Expression Host</b>	HEK293
<b>Tag</b>	N-6×His Tag
<b>Molecular Characterization</b>	6×His tag Mouse TENM4(Gly367-Arg2771)
<b>Molecular Weight</b>	The protein has a predicted molecular mass of 269.0 kDa after removal of the signal peptide.
<b>Purity</b>	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Involved in neural development, regulating the establishment of proper connectivity within the nervous system. Plays a role in the establishment of the anterior-posterior axis during gastrulation. Regulates the differentiation and cellular process formation of oligodendrocytes and myelination of small-diameter axons in the central nervous system (CNS). Promotes activation of focal adhesion kinase. May function as a cellular signal transducer.[UniProtKB/Swiss-Prot Function]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated






Figure 1. Mouse TENM4 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

