

PRODUCT INFORMATION

Uniprot ID	Q13641
Common Name	H8, mAb5T4
Conjugate	PE-conjugated
Synonyms	M6P1;WAIF1;TPBG
Applications	Flow Cyt
Recommended Dilutions	1µl/test
Formulation & Reconstitution	Liquid□PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	5T4
Description	PE-conjugated Anti-5T4 (H8 biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1 µl/test PE-conjugated Anti-5T4 (H8 biosimilar) mAb on HEK293 cells (Blue histogram) or Jurkat cells(Red histogram).

