

PRODUCT INFORMATION

Common Name	A06
Conjugate	PE-conjugated
Synonyms	PLAP-1;PLAP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - Lambda1
Reactivity	Human
Target	ALPP
Uniprot ID	P05187
Description	PE-conjugated Anti-ALPP(TCRCure A06 Biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

