

PRODUCT INFORMATION

Common Name	MG-A271, MGA271
Conjugate	PE-conjugated
Synonyms	B7-H3;CD276;B7 homolog 3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	B7-H3
Uniprot ID	Q5ZPR3
Description	PE-conjugated Anti-B7-H3 (enoblituzumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

