

PRODUCT INFORMATION

Common Name	4F12
Conjugate	PE-conjugated
Synonyms	PD1, PDCD1, CD279, SLEB2, hPD-1, hPD-I, hSLE1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Rabbit
IgG type	Rabbit IgG
Reactivity	Canine
Target	PD-1
Uniprot ID	NP_001301026.1
Description	PE-conjugated Anti-Canine PD-1(INTERVET 4F12) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

