

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | PRO 140, HuPRO-140, PRO-140 |
| Conjugate | PE-conjugated |
| Synonyms | CKR5; CCR-5; CD195; CKR-5; CCCKR5; CMKBR5; IDDM22; CC-CKR-5 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Humanized |
| IgG type | Human IgG4(S228P) - Kappa |
| Reactivity | Human |
| Target | CCR5 |
| Uniprot ID | P51681 |
| Description | PE-conjugated Anti-CCR5(leronlimab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

