

PRODUCT INFORMATION

Uniprot ID	P13726
Common Name	HuMax-TF
Conjugate	PE-conjugated
Synonyms	Tissue factor, TF, F3, Thromboplastin, Coagulation factor III
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG1 – kappa
Reactivity	Human
Target	CD142
Description	PE-conjugated Anti-CD142(tisotumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

