

PRODUCT INFORMATION

Common Name	FMC63
Conjugate	PE-conjugated
Synonyms	CD19;4;VID3;GC12802
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Mus musculus
IgG type	Mouse IgG2a - Kappa
Reactivity	Human
Target	CD19
Uniprot ID	P15391
Description	PE-conjugated Anti-CD19(FMC63 biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

