Cat. No. DMC101114P

PRODUCT INFORMATION

Clone ID 1D4 CD205 **Target**

Synonyms LY75;CLEC13B;DEC-205;GP200-MR6;LY-75

Host Species

PE-conjugated Anti-CD205 antibody(1D4), IgG1 **Description**

Chimeric mAb

Delivery Under Development

Uniprot ID 060449

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

Flow Cyt 1:100 **Dilutions**

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Reconstitution

Background

Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

Acts as an endocytic receptor to direct captured

antigens from the extracellular space to a specialized antigen-processing compartment (By

similarity). Causes reduced proliferation of B-lymphocytes.[UniProtKB/Swiss-Prot Function]

Research use only **Usage**

Conjugate PE-conjugated

All DIMA recombinant antibodies are genuinely

Email: info@dimabio.com Website: www.dimabio.com

generated by DIMA Biotech. They are all under patent application. Any protein sequencing or **DIMA Disclaimer** reverse engineering attempt is prohibited. We are

actively scrutinizing all patent application to ensure no IP infringement.

