

PRODUCT INFORMATION

Uniprot ID	P25942
Common Name	CFZ-533, NVP-CFZ533, OM11-62MF
Conjugate	PE-conjugated
Synonyms	CD40;Bp50;CDW40;MGC9013;TNFRSF5;p50
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	CD40
Description	PE-conjugated Anti-CD40 (iscalimab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

