

PRODUCT INFORMATION

Common Name	MEDI-9447, MEDI9447
Conjugate	PE-conjugated
Synonyms	NT5E, 5'-Nucleotidase, 5'-NT, NT5, NTE
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG1(K214R,L234F,L235E,P331S)
Reactivity	Human
Target	CD73
Uniprot ID	P21589
Description	PE-conjugated Anti-CD73(oleclumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

