

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | RO5541077-IG-7596 |
| Conjugate | PE-conjugated |
| Synonyms | AGM6;B29;IGB |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid-PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Chimeric Rabbit |
| IgG type | Rabbit IgG |
| Reactivity | Human |
| Target | CD79B |
| Uniprot ID | P40259 |
| Description | PE-conjugated Anti-CD79B (polatuzumab biosimilar) antibody, Rabbit IgG Chimeric mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

