

PRODUCT INFORMATION

Uniprot ID	P40259
Common Name	RO5541077-IG-7596
Conjugate	PE-conjugated
Synonyms	AGM6;B29;IGB
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid-PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Rabbit
IgG type	IgG
Reactivity	Human
Target	CD79B
Description	PE-conjugated Anti-CD79B (polatuzumab biosimilar) antibody, Rabbit IgG Chimeric mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

