

PRODUCT INFORMATION

Uniprot ID	P08195
Common Name	HuM23
Conjugate	PE-conjugated
Synonyms	SLC3A2; 4F2; MDU1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(K97R)
Reactivity	Human
Target	CD98
Description	PE-conjugated Anti-CD98(HuM23 biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

