

PRODUCT INFORMATION

Common Name	SAR-408377
Conjugate	PE-conjugated
Synonyms	CEA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	Human IgG1(E356D,M358L) - kappa
Reactivity	Human
Target	CEACAM5
Uniprot ID	P06731
Description	PE-conjugated Anti-CEACAM5(tusamitamab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

