

PRODUCT INFORMATION

Common Name	AMG 334,AMG-334,erenumab-aooe
Conjugate	PE-conjugated
Synonyms	CALCRL
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:400
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG2 - Kappa
Reactivity	Human
Target	CGRPR
Uniprot ID	Q16602
Description	PE-conjugated Anti-CGRPR(Erenumab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 100 μ l/test (1:400) PE-conjugated Anti-CGRPR(Erenumab biosimilar) mAb (BME100166P) on HEK293 cells transfected with Human CGRPR and RAMP1 protein (Blue histogram) and HEK293 cells (Red histogram) .

