Cat. No. DMC100673P

## **PRODUCT INFORMATION**

**Clone ID DMC673 CRTAM Target** CD355 **Synonyms Host Species** Rabbit

PE-conjugated Anti-CRTAM antibody(DMC673); **Description** 

IgG1 Chimeric mAb

**Delivery Under Development** 

**Uniprot ID** 095727-1

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

**Background** 

Flow Cyt 1:100 **Dilutions** 

Purified from cell culture supernatant by affinity **Purification** 

chromatography

Formulation & Reconstitution

Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

The CRTAM gene is upregulated in CD4 (see MIM 186940)-positive and CD8 (see CD8A; MIM 186910)-positive T cells and encodes a type I transmembrane protein with V and C1-like Ig domains (Yeh et al.; 2008 [PubMed 182203701) [supplied by CMIM: Feb 2000]

18329370]).[supplied by OMIM; Feb 2009]

**Usage** Research use only

Conjugate PE-conjugated

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or

> > Email: info@dimabio.com Website: www.dimabio.com

**DIMA Disclaimer** reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to

ensure no IP infringement.

