

PRODUCT INFORMATION

Common Name	DKN-01, LY-2812176
Conjugate	PE-conjugated
Synonyms	DKK-1;SK
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	DKK1
Uniprot ID	O94907
Description	PE-conjugated Anti-DKK1(sirexatamab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

