

PRODUCT INFORMATION

Uniprot ID	P00533
Common Name	QD6-AZD9592 without cMet mAb
Conjugate	PE-conjugated
Synonyms	ERBB; ERBB1; HER1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid-PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	EGFR
Description	PE-conjugated Anti-EGFR(QD6 biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

