

PRODUCT INFORMATION

Common Name	XENP28278
Conjugate	PE-conjugated
Synonyms	CD203c;NPP3;PD-Ibeta
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:400
Formulation & Reconstitution	Liquid \square PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	ENPP3
Uniprot ID	O14638
Description	PE-conjugated Anti-ENPP3(Xencor) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMA BIO CONFIDENTIAL






Figure 1. Flow cytometry analysis with 100 μ l/test (1:400) PE-conjugated Anti-ENPP3(XENP28278) mAb (BME100151P) on HEK293 cells transfected with human ENPP3 (Blue histogram) or HEK293 cells (Red histogram).

DIMA BIO CONFIDENTIAL

