

PRODUCT INFORMATION

Uniprot ID	O14638
Common Name	XENP28278
Conjugate	PE-conjugated
Synonyms	CD203c;NPP3;PD-Ibeta
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:400
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	ENPP3
Description	PE-conjugated Anti-ENPP3(Xencor) mAb
Delivery	In Stock
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL






Figure 1. Flow cytometry analysis with 100 μ l/test (1:400) PE-conjugated Anti-ENPP3(XENP28278) mAb (BME100151P) on HEK293 cells transfected with human ENPP3 (Blue histogram) or HEK293 cells (Red histogram).

DIMA BIO CONFIDENTIAL

