

PRODUCT INFORMATION

Uniprot ID	Q96RD9
Common Name	BFCR-4350-A, BFCR-4350A, RG-6160, RO-7187797
Conjugate	PE-conjugated
Synonyms	BXMAS1;CD307;CD307e;FCRH5;IRTA2;PRO820
Applications	Flow Cyt
Delivery	In Stock
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	FCRL5
Description	PE-conjugated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb
Recommended Dilutions	Flow Cyt 1:200
Storage & Shipping	Store at 2°C-8°C for 6 months
Formulation & Reconstitution	Liquid□PBS with 0.1% Proclin300, 0.18% BSA
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scr






Figure 1. Flow cytometry analysis with 100 μ l/test (1:200) PE-conjugated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb (BME100089P) on HEK293 cells transfected with human FCRL5 (Blue histogram) or HEK293 cells(Red histogram).

