

PRODUCT INFORMATION

Uniprot ID	P47871
Common Name	AMG 477, REMD-477
Conjugate	PE-conjugated
Synonyms	Glucagon receptor;GL-R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG2
Reactivity	Human
Target	GCGR
Description	PE-conjugated Anti-GCGR(volagidemab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

