

PRODUCT INFORMATION

Uniprot ID	Q99988
Common Name	PF-06946860
Conjugate	PE-conjugated
Synonyms	GDF-15;MIC-1;MIC1;NAG-1;PDF;PLAB;PTGFB
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	GDF15
Description	PE-conjugated Anti-GDF15(ponsegromab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scr

