

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | CAM-3001 |
| Conjugate | PE-conjugated |
| Synonyms | GM-CSF-R-alpha; GMCSFR-alpha; GMR-alpha; CDw116; CD116; CSF2RA; CSF2R; CSF2RY |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Homo sapiens |
| IgG type | Human IgG4(S228P) - Kappa |
| Reactivity | Human |
| Target | GM-CSFR |
| Uniprot ID | P15509 |
| Description | PE-conjugated Anti-GM-CSFR(mavrilimumab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

