

PRODUCT INFORMATION

Uniprot ID	Q9NZD1 and P07766
Common Name	LBL-034
Conjugate	PE-conjugated
Synonyms	GPRC5D; CD3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	GPRC5D and CD3
Description	PE-conjugated Anti-GPRC5D and CD3(LBL-034 biosimilar) BsAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

