

PRODUCT INFORMATION

Common Name	LBL-034
Conjugate	PE-conjugated
Synonyms	GPRC5D
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1(L234A,L235A,D265A,P329A,Y349C,E356D,M358L,T366S,L368A,Y407V) - Kappa
Reactivity	Human
Target	GPRC5D
Uniprot ID	Q9NZD1
Description	PE-conjugated Anti-GPRC5D(LBL-034 biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

