

PRODUCT INFORMATION

Uniprot ID	Q9NZD1
Common Name	LBL-034
Conjugate	PE-conjugated
Synonyms	GPRC5D
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(L117A, L118A, D148A, P212A, Y232C, E239D, M241L, T249S, L251A, Y290V )
Reactivity	Human
Target	GPRC5D
Description	PE-conjugated Anti-GPRC5D(LBL-034 biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

