

PRODUCT INFORMATION

Common Name	5F9,MLN2045
Conjugate	PE-conjugated
Synonyms	DIAR6;GC-C;GUC2C;MECIL;MUCIL;STAR
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1µl/test
Formulation & Reconstitution	Liquid□PBS with,0.18%BSA□0.1%Proclin300
Host Species	Homo sapiens
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	GUCY2C
Uniprot ID	P25092
Description	PE-conjugated Anti-GUCY2C(indusatumab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1 μ l/test PE-conjugated Anti-GUCY2C(indusatumab biosimilar) mAb on HT55 cells (Blue histogram) or CHO-S cells (Red histogram).

